

Esophagus Bougie Dilator



Diameter: ϕ 5, 6, 7, 8, 9, 10, 11, 12,13, 14,15, 16, 17, 18, 19, 20.


Length: 700mm or 850mm.

Package: Single-use dilators choose tyvek pouch package; reusable sets choose alloy box package.

Features:

- It is used for the expansion and treatment of esophagus and cardia.
- Made by medical PVC materials with proper softness.
- The front cone angle is manufactured by mechanized processing. With smooth transition, it can be inserted into the narrow parts successfully.

Disposable Biopsy Forceps



REF	Jaw Dia. (mm)	Working Length (mm)
HQ-1.8x1050-*-*	1.8	1050
HQ-1.8x1100-*-*	1.8	1100
HQ-1.8x1200-*-*	1.8	1200
HQ-1.8x1600-*-*	1.8	1600
HQ-1.8x1800-*-*	1.8	1800
HQ-2.3x1600-*-*	2.3	1600
HQ-2.3x1800-*-*	2.3	1800
HQ-2.3x2300-*-*	2.3	2300

Note: please indicate " with needle (W)", "no needle (O)" , "plastic coated (P)", "no plastic coated (O)" .

Disposable Endoscopic Ligating Loop



REF	Suture Diameter (mm)	Length of Plastic Tube (mm)
TZ-300	0.4-0.5	305

Features:

- Mostly used in laparoscopic ligation surgery, gall bladder surgery, bile duct surgery, fibroid surgery, fallopian tube surgery, appendix surgery, intestine surgery etc.
- Polypropylene suture (non-absorbable) has higher intensity and makes stronger knot.