


Disposable Hemorrhoid Multi-band Ligature with Anoscope



Applicable Operation

Applicable for internal hemorrhoid in all stages, internal hemorrhoid of mixed hemorrhoid, anal papillae and rectal polyp.



Production Characteristics

- 1 Realized the automation, the ligation treatment process to save time, effort, practical, and easy;
- 2 Patient safety, avoid to take medicine, no injection, no cross infection;
- 3 Single complete diagnosis and surgical operation, only take 5 to 10 minutes;
- 4 After treatment, the complications are rare.

Technological Parameters

Model	Specification	Colonoscopy Specification	Outer diam Adsorption tube	Inner Diam Ligation ring	Outer Diam Ligation ring	Outer diam Colonoscopy	Height of Ligation ring
RYTZ	4	-(A/B)	12	2.0	5.0	20	1.8
	5						
	6						
	7						